




BUILD YOUR HYBRID

HEALTHCARE INTEGRATION



# MEDICAL APPLICATION CONNECTION HUB – MACH

MACH aims to address the integration and data challenges faced by the healthcare industry by providing hybrid cloud integration for health data. MACH is a secure, open standards-based, cloud-hosted platform for quick and simple deployment of on-premises and cloud-native integration services. MACH allows you to manage instance creation, updates and replacement centrally. It is highly scalable, API-driven, and privacy-first by default. MACH seeks to

address every use case and serve as a one-stop shop for partners facing integration challenges. MACH will be made available as a SaaS model to cover the entire lifecycle of an interface. Alternatively, it can also be acquired as an intelligent DaaS architecture solution that can be controlled via MACH to securely transfer data even in large amounts or with a high-velocity and at various QoS levels.


MACH provides a secure, reliable and compliant platform for sharing health data between medical applications.

The Orchestra Health Service Bus (OHSB) on which MACH is built provides industry standards-based message transformation tools for all flavors of HL7 (v2 and v3), FHIR<sup>1</sup>, and supporting standards such as ASTM, XDT, EDIFACT, CDA, CCD, CCDA and DICOM. MACH also supports device integrations<sup>2</sup> and all widely used communication technologies, such as MLLP, TCP/IP, JDBC, SOAP, REST, DICOM Worklist, and TLS certificate-based exchange. The built-in DICOM routing features offer DICOMweb and DIMSE, MPPS, plus change and (semi-automated) anonymization functions, which can be implemented as pre- and post-processing steps facilitated by MACH. MACH also offers a hybrid configuration: if required, you can customize the centrally managed cloud platform by connecting to an on-premise instance.

How you manage and monitor data flows determines how effective a data integration platform is. On the MACH platform, monitoring is possible from a cluster perspective all the way down to individual instances, allowing for high-level status investigation as well as in-depth scenario exploration and log search capabilities using tools such as Grafana.

The use case illustrated next uses MACH's Cloud Connector (CC) to transfer DICOM studies from an on-premises image source to a cloud-hosted AI solution. The fully fledged Health Service Bus that is included ensures that MACH CC functions flawlessly and can securely and reliably transfer protected health information (PHI<sup>3</sup>) right out of the box, with optional pre- and post-processing steps.

Fast, frictionless and secure transfer of big data through data streaming in MACH Cloud Connector (CC)



MACH can be used in a variety of ways. Due to the modularity of the various options, customers benefit from additional flexibility and

can tailor the architecture to their individual needs.

Standard MACH deployment options include:<sup>4</sup>

- **Lightweight**  
 Docker or VM based image delivery, local database therefore limited data resilience suitable for on-premise and cloud use. Will be managed from a central point and receive pushed configuration changes.
- **Default**  
 Micro-Kubernetes Cluster with shared databases, including cloud backup suitable for on-premise and cloud use. The default option can be provided in multiple variations depending on the focus or number of scenarios. It can be understood as a mid- to high-range interoperability layer capable of handling heavy loads.
- **High availability and cluster**  
 Zero downtime and high availability is the main goal of implementations in a k8s cluster providing a robust operations-centric architecture that is highly scalable and resilient. Offering container self-healing capabilities and designed to efficiently manage large amounts of containers that are distributed across complex cloud environments.

Basic deployment options with close-to-zero downtime are available for all MACH instances. Rolling deployments can be provided in a container environment such as Docker or k8s.

Additional deployment scenarios can be implemented based on the client's needs.

<sup>1</sup> HL7, FHIR and the  Logo are registered trademarks belonging to Health Level Seven International

<sup>2</sup> Native device integration requires a hardware dongle from the vendor to transform serial inbound stream

<sup>3</sup> Protected Health Information (HIPAA Administrative Simplification Provision)

<sup>4</sup> There are prerequisites for on-premise deployment: either docker runtime or VM Hypervisor on-premise deployment needs to be available at the customers site

# Share your pressing integration challenges with us!

Get in touch.

**Our hybrid integration experts are happy to help:**

+44 203 983 9860 | [info@mach.health](mailto:info@mach.health) | [mach.health](https://mach.health)

**x-tention Limited**

Bourne Gardens, Exeter Park Road  
Bournemouth, Dorset  
United Kingdom BH2 5BD

**xtention**  
IT with care.